

Medical Measuring Systems and Scales for

# PHYSICIANS



## seca 216

Mechanical measuring rod for children and adults

- + Easy-to-read measurement scale
- + Adjustable head positioner
- + Choice of two units of measurement
- + Simple assembly



Measuring range	max. 138 cm / 55.3" selectable from 3.5 – 230 cm / 1.4 – 90"
Graduation	1 mm / 1/16"

## seca 206

Mechanical measuring tape

- + Space-saving, roll-up mechanism
- + Simple read-out of results through display window
- + Durable tape made of metal



Measuring range	0 – 220 cm or 0 – 87"
Graduation	1 mm or 1/8"

You can feel seca quality, it gives a good impression and pays off. The caliper gauge of a seca stadiometer runs smooth and with millimetre precision over the abrasion-resistant scale and a seca column scale gives a safe and high-quality impression, even after thousands of weighings. It is not without reason that seca is the world market leader in the field of medical weighing and measuring and thus the quality standard for every medical practice.



Image contains optional digital measuring rod seca 234 and barcode scanner from a third-party vendor.



### seca 354 m / 354

#### 2-in-1 mobile baby scale and flat scale for toddlers

- + Weigh babies lying down and toddlers standing up
- + Easily detach the weighing tray from the base to use as a flat scale
- + High-quality baby scale with calibration approval at an economical price
- + Wide range of functions: TARE, HOLD and Breast-Milk-Intake-Function (BMIF)



Capacity	20 kg / 44 lbs
Graduation	5 g < 10 kg > 10 g / 0.2 oz < 22 lbs > 0.5 oz
Power supply	seca 354 m: Power adapter, batteries seca 354: Batteries
Functions	TARE, BMIF, Auto-HOLD, kg/lbs switch-over, automatic switch-off
Accessories	+ Mobile carts seca 403 and seca 402 + Measuring board seca 417 or measuring mat seca 210 + Backpack seca 409 (for combination with measuring board or measuring mat) + Carrying case seca 414 (for combination with measuring board) + Backpack seca 431 + Carrying case seca 413

### seca 333 i

#### EMR-validated baby scale with Wi-Fi function

- + Damping system for fast weighing
- + Easy-to-read, backlit LCD display
- + USB port for barcode scanner for user and patient identification




Capacity	20 kg / 44 lbs
Graduation	5 g / 0.2 oz
Power supply	Power adapter, batteries
Functions	TARE, BMIF, HOLD, Auto-HOLD, damping, automatic switch-off, Wi-Fi, kg/lbs switch-over, confirm button
Ports	USB for barcode scanner
Accessories	+ Mobile carts seca 403 and seca 402 + Digital measuring rod seca 234 + Analog measuring rod seca 232 n + Battery for longer mobile use seca 454 + Carrying case seca 428

Image contains optional handrail seca 477.

### seca 777

#### Digital column scale with eye-level display

- + Medical-grade, robust scale designed for years of use with high patient volume
- + Built-in measuring rod with 10 – 230 cm measuring range
- + USB interface for EMR integration solutions or PC connectivity
- + Large rubber coated wheels for excellent movability
- + Optional: Handrail seca 477 provides extra patient security
- + Optional: Panda bear facilitates measuring and weighing of children



<b>Scale</b>	
Capacity	250 kg / 550 lbs
Graduation	100 g / 0.2 lb
Functions	TARE, HOLD, kg/lbs switch-over lock, BMI
Accessories	+ Handrail seca 477 + Panda seca 459 + Power adapter 401

<b>Measuring rod (included)</b>	
Measuring range	10 – 230 cm / 4 – 90"
Graduation	1 mm / 1/8"



### seca 874 dr

#### Flat scale with customizable label

- + Customizable label helps minimize scale misplacement
- + Inexpensive, reliable and durable: seca quality that pays off
- + Double display allows the doctor and patient to read simultaneously
- + Mother+child function for hassle-free weighing of toddlers

[seca.com/doctorscale](http://seca.com/doctorscale)



Capacity	200 kg / 440 lbs
Graduation	50 g < 150 kg > 100 g 0.1 lb < 330 lbs > 0.2 lb
Functions	Mother+child function, Auto-HOLD, kg/lbs switch-over, automatic switch-off
Accessories	+ Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element) + Power adapter 447 + Switch-mode power adapter 400